Hall Effect

Lab 5

Dhruva Sambrani, Gowri A (Group 6)

April 11, 2022

Contents

Aim	:	
Method	2	
Observation	2	
Results	3	

Aim

1. To measure the Hall effect in a semiconductor, and hence derive certain electrical properties of the semiconductor.

Method

- 1. Measure the thickness of the sample as t and resistance as R
- 2. Make connections such that you can pass current and measure the voltage perpendicular to each other in the sample.
- 3. Connect and place the Electromagnet such that the resulting magnetic field is perpendicular to both current and voltage axes.
- 4. Place a gauss meter between Electromagnet and sample
- 5. Calibrate Gauss meter
- 6. Pass current through the Electromagnet such that the gauss meter reads ~ 1.5 k gauss. Record it as B
- 7. Pass current through the sample from about 1mA to 15 mA.
- 8. Measure the corresponding voltage.
- 9. Plot and fit the graph to y = mx, with m giving us mean value of V_H/I
- 10. Use m to find the following-
 - 1. $R_H = mt/B$
 - 2. $n = 1/R_H e$
 - 3. r = Rt
 - 4. $\mu = R_H/r$
 - 5. $\theta_H = \tan^-(\mu B)$

Observation

Quantity	Value
t	5e-4m
В	0.1464 T
\mathbf{R}	0.870Ω

V (mV)	I (mA)	V (mV)	I (mA)
0.86	8.1	16.7	165.5
1.93	18.3	17.14	171.2
3.02	28.5	-1.91	-17.8
4.06	38.2	-3.42	-32.1
5	47.1	-4.46	-42.1
6	56.9	-5.16	-48.6
7.35	69.6	-7.06	-67.1
7.99	75.8	-8.55	-81.6
9.42	89.8	-9.28	-88.8
10.16	97.3	-10.45	-100.6
11.11	106.8	-11.84	-114.9
12	116	-12.97	-126.8
13.16	128	-15.56	-153.9
14.71	144.2	-17.07	-170.1
15.6	153.8		




Figure 1: Hall Effect Experimental Setup

Results




Figure 2: Plot

- 1. Mean of $V_H/I=9.76\pm0.3$
- 2. Hall Coefficient $R_H = 0.0333 \pm 2.368$
- 3. $n = 1.88 \cdot 10^20 \pm 2.368$
- 4. Resistivity $r = 4.35 \cdot 10^-4 \pm 2.11$
- 5. Mobility $\mu = 76.6 \pm 4.478$
- 6. Hall angle $\theta_H = 1.48 rad$